

Successes and Challenges in Applications of a Laboratory-Based Scanning XPS/HAXPES Instrument

K. Artyushkova*¹, A. Vanleenhove², T. Conard², P. -M. Deleuze³ and O. Renault³,

¹Physical Electronics Inc., 18725 Lake Drive E, Chanhassen, MN 55317, USA


²IMEC, 3001 Leuven, Belgium

³Univ. Grenoble-Alpes, CEA, Leti, 38000 Grenoble, France



Introduction

With the recent development of lab-based Hard X-ray Photoelectron Spectrometers, there are exciting new opportunities for looking beneath the surface layers accessible using conventional soft X-ray sources in XPS. We describe the development and applications of a laboratory-based instrument, the PHI Quantes, equipped with two scanning microprobe monochromated X-ray sources, Al K α (1486.6 eV) and Cr K α (5414.8 eV), thus enabling both traditional XPS and HAXPES experiments in the same instrument.



Left, schematic of X-ray source and analyzer on PHI Quantes. Software-controlled anode and shutter can automatically switch between Al and Cr X-rays in < 1 minute. Right, an overview photo of PHI Quantes.

Benefits of HAXPES

- Deeper depth of analysis:**
 - Buried interfaces:
 - many interfaces such as thin films on a substrate
 - active electronic layers below a surface capping layer
 - depth-profiling through heterostructures and e.g. layered low-dimension materials
 - probing of dopants and contaminants in the bulk of a material
 - Reduced effects of surface contamination
 - Reduced effects of ion beam damage at the interface in depth profiles
- Access to higher energy core photoelectron peaks:**
 - Additional transitions of varied kinetic energy and hence sampling depth
- Shifted Auger Transition:**
 - Avoid overlap between core electron lines and Auger peaks
 - Use of Auger Parameter

Analysis of Al₂O₃/GaN thick gate oxides

Challenge: Analyzing thick oxide layers


Benefits of using HAXPES:

- Probing beneath the thick oxide layer without sputtering
- Avoiding critical overlap between N 1s and Ga LMM Auger line when using Al K α source.

Goal of analysis:

- Probe GaN layer beneath a 30 nm oxide layer which is identical to that used in MOSc-HEMTs devices
- Analyze samples with increasing thickness from 5 till 30 nm


Comparison of the Cr K α HAXPES Ga2p_{3/2} peak intensity for different Al₂O₃ thicknesses.



Thickness, nm

nominal	HAXPES
5 nm	4.5
10 nm	8.2
20 nm	18.9
30 nm	27.8

Structure of 20 nm Al₂O₃ film modeled using StrataPHI



HC: 1.4 ± 0.6 nm
Al₂O₃: 18.9 ± 1.0 nm

Success:

- Ga was detected even for 30 nm oxide layer sample
- The thickness of alumina layers estimated using StrataPHI software for thin-film structure optimization is within 10% of the nominal value.

Damage-free depth profiling

Challenge: Probing beneath ion beam damage

Sampling depth using different X-ray sources, nm (highlighted depth > 15 nm)

- For efficient and accurate chemical analysis of buried layers and interfaces, the photoelectrons analyzed must originate from below the ion beam damage depth, from depths of ~15 nm or more.




BE, eV	Al	Ag	Cr
Ga 2p	1120	3.4	11.1
Zn 2p	1022	3.0	8.7
Ni 2p	850	3.3	8.3
Sn 3d	486	7.6	15.5
Ti 2p	458	6.6	13.3
C 1s	285	9.2	17.6
Si 2p	100	9.2	16.6
Al 2p	74	9.3	16.7
Hf 4f	17	8.4	14.4

Benefits of using HAXPES:

- Probing beneath the ion beam damage depth

Goal of analysis:

- For a direct comparison of the decreased sensitivity to ion beam damage using a Cr source versus an Al source, acquire depth profile from a stack of 30 nm TiN/ ~5 nm SnO₂ films on a SiO₂/Si substrate using both X-ray sources.



The spectra shown are at 30, 31, and 37 minutes for the Al profile and at 21, 22, and 33 minutes for the Cr profile

Success:

- The degree of reduction observed in HAXPES data at the maximum concentration of Sn is much smaller than in XPS data obtained using the Al X-ray source.

Experimental elemental RSF database


Challenge: Accurate quantification of HAXPES data

Protocol:


- RSF of transitions was derived based on the ratio of intensities with respect to pure silver

$$RSF_{transition} = \frac{Intensity_{transition}}{RSF_{Ag3d}} = \frac{TFC_{transition} \cdot Asymmetry_{transition}}{TFC_{Ag3d} \cdot Asymmetry_{Ag3d}}$$

- RSF_{Ag3d} wrt to RSF_{F1s} = 1 was derived by analyzing Ag triflate and ionic liquid on three different instruments



Cr K α reference spectra



- Reference compounds providing necessary chemical shift information for newly observed transitions

Beta database v2 released:

- Transitions for 52 elements
- Reference ionic liquid for F, C, O, N and S
- Pure bulk metals for all others

Success:

- Accurate quantification using all transitions

	Nominal	Si1s	Si2s	Si2p
SiO ₂	33.3%	33.1%	31.8%	32.9%
SiC	50.0%	49.4%	47.9%	50.9%
Si ₃ N ₄	42.9%	42.8%	39.6%	42.5%
		Al1s	Al2s	Al2p
Al ₂ O ₃	40.0%	38.5%	36.2%	37.1%
AlN	50.0%	51.4%	50.0%	51.2%
		Mo/Se	2s/2s	2p/2p
MoSe ₂	33.3%	35.7%	33.8%	31.6%

Analysis of 2D materials


Challenge: Effect of surface contamination on quantification

Benefits of using HAXPES:

- Probing beneath the adventitious carbon and oxygen without using ion beam sputtering

Goal of analysis:

- M_{n+1}AX_n (MAX) are precursors for MXenes – 2D transition metal carbides and nitrides
- Use Cr X-ray microprobe to analyze individual MAX particles of proper orientation
- Quantify oxygen defects using Cr source



5x5 mm² mosaic of stitched X-ray induced secondary electron images (SXI)

Quantification results

theory	Ti/Al=3	Ti/C=1.5	O1s
Al (4-6 nm)	0.70	0.27	51.6
Cr (16-20 nm)	3.28	1.70	29.6

Success:

- As-expected stoichiometry is seen from HAXPES data
- A significant amount of oxygen is present at deeper sampling depths, suggesting oxygen defects within the MAX structure