



Aroclors, Homologs and Congeners – An Evaluation of the Options for PCB Analysis and a Comparison of the Interpretive Value

Charles Neslund, Scientific Officer

TCEQ Environmental Trade Fair, May 15-16, 2018

*Providing comprehensive scientific resources to
environmental clients worldwide.*

Polychlorinated Biphenyls - PCBs



- Most commonly known for use in transformers as dielectric fluids
- Also used in hydraulic fluids, lubricants, as plasticizers and in caulking compounds
- Primary manufacturer was Monsanto – product given the trade name Aroclor
- Aroclor nomenclature indicated the average % chlorine in the mixture
i.e. Aroclor 1254 = approx 54% chlorine

Why are they important?



- PCBs bioaccumulate
- WHO (World Health Organization) has determined that some of the 209, ones referred to as “dioxin-like” PCBs, dl-PCBs, or coplanar PCBs, may be anticipated to cause cancer
- People can be exposed to PCBs through breathing in contaminated air, consuming contaminated food, and by skin contact with old electrical equipment that contains PCBs.

Polychlorinated Biphenyls - PCBs



PCBs can be grouped or identified by three common descriptions;


- Aroclor – a mixture of up to one hundred different congeners, grouped together based on an average percent weight chlorine basis
- Homolog – congeners of the same chlorination level. There are 10 chlorination groupings
- Congener – any of the individual chlorinated biphenyl compounds of which 209 are possible

Polychlorinated Biphenyls - Aroclors




- The most common analysis technique is GC/ECD, SW-846 Method 8082
- Identifies and quantifies PCBs as an Aroclor
 - Aroclor 1016 Aroclor 1221 Aroclor 1232
 - Aroclor 1242 Aroclor 1248 Aroclor 1254
 - Aroclor 1260 Aroclor 1262
- Limits are typically in the ug/kg and ug/l range
- Technique is prone to interferences


Polychlorinated Biphenyls - Aroclors



Aroclor 1242



Aroclor 1254



Aroclor 1260

Polychlorinated Biphenyls - Aroclors



- Each Aroclor mix consists of dozens of congeners
- Congeners with the same number of chlorines on the biphenyl ring are called homologs
- For example, Aroclor1254 contains;
 - 1% Trichlorobiphenyls
 - 15% Tetrachlorobiphenyls
 - 53% Pentachlorobiphenyls
 - 26% Hexachlorobiphenyls
 - 4% Heptachlorobiphenyls

Why are they important?



- Everything you know about Aroclor PCBs applies in some part to PCB Congeners
- Not all PCB congeners are found in Aroclors
- Approximately 130 of the 209 congeners are detected in Aroclors
- Congeners are also categorized by the number of chlorines substituted on the biphenyl ring
- These categories are called homolog groups or isomer groups

Polychlorinated Biphenyls - Homologs



Homologs can be characterized by EPA Method 680 or EPA Method 1668.

EPA Method 680

- Gas Chromatography/Mass Spectrometry (low resolution)
- Can be operated full scan or selected ion monitoring (SIM)
- Therefore, sensitivities similar to what you have with 8270

0.1 ug/l for the lower chlorination homologs

0.3-0.5 ug/l for the higher chlorination homologs

Low/single digit ug/kg for soils/solids

Polychlorinated Biphenyls - Homologs




Nine congeners are used to “characterize”
the homolog groups

PCB 1	Monochlorobiphenyls
PCB 5	Dichlorobiphenyls
PCB 29	Trichlorobiphenyls
PCB 50	Tetrachlorobiphenyls
PCB 87	Pentachlorobiphenyls
PCB 154	Hexachlorobiphenyls
PCB 188	Heptachlorobiphenyls
PCB 200	Octachlorobiphenyls
PCB 209	Nona- and decachlorobiphenyl

Total Ion Chromatogram of 680 CCAL



Abundance



Time->

Polychlorinated biphenyls - Homologs




- Quantitation is achieved with internal standard calibration using deuterated PAHs, like 8270
- Homolog group is determined by summing the areas of all peaks within a given retention time, that have the same mass
- Even though a select few congeners are used for the calibration response, data does not yield congener information
- Results will consist of a value for each of the homolog groups and a sum of the homolog groups for a PCB total



TIC of Typical Sample

Abundance




Time→

TIC of Contaminated Sample




Abundance




Extracted Ion for Pentachlorobiphenyl



Abundance



Time →
Abundance




Time →

What are PCB Congeners?



A congener is a specific chlorine substitution pattern on a biphenyl ring.

There are a total of 209 possible substitution patterns ranging from 1 chlorine up to 10 chlorines



Polychlorinated Biphenyls – Congeners



- To monitor and analyze for the more toxic congeners, High Resolution Mass Spectrometry (HRMS) is employed
- EPA Method 1668
 - HRGC/HRMS
 - Capillary column to resolve many of the congeners
 - Special clean-up techniques
- Reporting limits of ng/kg and pg/l

Polychlorinated Biphenyls – Congeners



- Able to calibrate for all 209 congeners
- Calibration is also mostly internal standard except that the compounds used as internal standards are isotopically labeled PCB congeners
- Able to differentiate between congeners that are constituents of Aroclors and those that are not.
- Able to differentiate contribution of dl-PCBs from total PCB content


Calibration Standard with 209 Congeners



Z:\16JUL13\16JUL13-02

7/13/2016 10:53:04 AM VER-CAL20941637B
SPB-Octyl 30 M x 0.25um x 0.25mm

RT: 9.00 - 63.00



NL:
1.17E6
TIC MS
16JUL13-
02

Example Data – PCB15




16JUL06-02 printed 7/6/2016 13:46
Sample VER-CAL-209416378 / PCBs209-0
Inst ID: DF19780-16JUL06 / Client:



Lancaster
Laboratories

Chromatogram

RT: 13.61 - 23.61 SM: 5G



Entry: N2PCB15 IS: L2PCB15



Lancaster Laboratories
Environmental

Polychlorinated Biphenyls – Congeners



- Can do fingerprinting of PCB sources based on ratios and relative amounts of the individual congeners
- Comparison of Aroclor based PCB congeners versus presence of non-Aroclor based PCB congeners
- Homolog results determined by summing of individual congeners that make up each chlorination level.
- Total PCB concentration as well as special case congeners quantified.

Special Case Congeners



Dioxin-like PCBs (12)

- PCB77
- PCB81
- PCB105
- PCB114
- PCB118
- PCB123
- PCB126
- PCB156
- PCB157
- PCB167
- PCB169
- PCB189

Indicator PCBs (6)

- PCB28
- PCB52
- PCB101
- PCB138
- PCB153
- PCB180

Questions

