

Navigating Rapid AST: Key Considerations for Direct-from-Blood Culture Implementation

April Bobenchik, PhD, D(ABMM)

Penn State Hershey Medical Center

abobenchik@pennstatehealth.psu.edu



PennState Health



PennState
College of Medicine

Disclosures

- Research funding:
 - BD, BioMérieux, Q-linea



PennState Health




PennState
College of Medicine

Objectives

- Summarize manual and currently available commercially rapid, direct-from-blood AST devices
- Identify key points/ considerations for device/system selection
- Outline some considerations and common challenges with implementing rapid AST



Approach To Implementation



PennState Health



PennState
College of Medicine

Antimicrobial Stewardship

Meet with the team

- ID pharmacists
- ID physicians
- Microbiologists
- IT / LIS

Identify the need

- Escalation
- De-escalation
- IV to PO
- Specific population

Availability to respond

- Day shift
- Nights and Weekends

Who will respond

- Stewardship member
- Pharmacist
- Infectious Disease
- Resident/ Fellow
- Primary team



PennState Health



PennState
College of Medicine

Clinical Stakeholders

Primary Team

- Responsible for acting on the results and making a change

Infectious Disease

- Consulted for difficult cases and therapeutic management

Quality

- Team following hospital quality metrics
 - Length of stay
 - Hospital readmissions

Hospital Leadership

- Help facilitate /champion new initiatives



PennState Health



PennState
College of Medicine

Commercial Manufacture

Review Website

- FDA status
- Product information
- Antimicrobials/ Organisms
- Posters/Papers/ Presentations

Attend Vendor Show

- Size
- Configuration
- Demo

Contact Local Sales Rep

- Specific Information
- Package Insert/ IFU
 - Limitations
 - Confirmatory testing
- Verification plan



PennState Health



PennState
College of Medicine

Information Technology/ Laboratory Information System

Laboratory Information System (LIS)

- Connectivity/
interface between
instrument and LIS
- Test build

Electronic Medical Record (EMR)

- What the result report will
look like to the clinician
- Special alerts
- Results for selected users
 - Stewardship




PennState Health



PennState
College of Medicine

Laboratory Considerations



PennState Health



PennState
College of Medicine

Ability to Obtain a Rapid Identification

Molecular ID



Verigene



Biofire Filmarray



Cobas eplex

- Result <2 h
- Limited to on-panel targets

MALDI TOF MS



Bruker MALDI Biotyper



BioMerieux VITEK MS

- Direct from blood culture
- Early growth on agar




PennState Health



PennState
College of Medicine

FDA clearance for a subset of Enterobacterales and Non-fermenters

Most Frequent Organisms



Useable Results

- ~ 85% Gram-negative bacilli are on panel
- ~ 15% of GNB cannot be resulted on a commercial rapid AST system
 - Alternative rapid AST (ie Disk diffusion) for monomicrobial off-panel Enterobacterales
 - Conventional AST

Banerjee et al, CID 2021
Ostermann et al, JCM 2024
Tibetts et al, JCM 2022
Snyder et al, JCM 2024
Esse et al, JCM 2023




PennState Health



PennState
College of Medicine

Need for Monomicrobial Blood Cultures

- ~ 3% polymicrobial GNB
 - Cannot distinguish on Gram Stain
 - Not always identified by rapid ID methods
 - Discovered upon sub-culture



- What you can do in the lab
 - Check your sub plate before releasing rapid AST results
 - Could do an early read before the rAST results are released
 - Review AST results before releasing to the chart
 - Look for unusual resistance patterns

Banerjee et al, CID 2021
Ostermann et al, JCM 2024
Tibetts et al, JCM 2022
Esse et al, JCM 2023




PennState Health



PennState
College of Medicine

Workflow Considerations




PennState Health



PennState
College of Medicine

Workflow Considerations



PennState Health



PennState
College of Medicine

Workflow Considerations

- Where are blood cultures incubated
 - On-site
 - Centralized lab
- How is the laboratory staffed
 - All 3 shifts
 - 1st and 2nd
 - 1st only
- Who is reading GS, performing rapid ID and AST
 - What shift, dedicated microbiologist or generalist
- Who is reporting the rapid AST result
 - Dedicated microbiologists, level of training



Clinical Considerations



How will the clinician view the result

Hyperlink to the instrument-generated report

Result View



12/21/2022 09:15

Abnormal

Culture Blood First...

Scans on Order 512414395

> Lab Result Scan - Scan on 12/22/2022 2:07 PM

IDENTIFICATION RESULTS

SPECIES	IDENTIFICATION
Klebsiella spp.	Positive ^{1, 2}

SUSCEPTIBILITY RESULTS

Klebsiella spp.

ANTIMICROBIAL	MIC	SIR
Amikacin	≤4	S ^{2, 3}
Ampicillin-Sulbactam	4	S ^{2, 3}
Aztreonam	2	S ^{2, 3}
Cefepime	≤1	S ^{2, 3}
Ceftazidime	≤2	S ^{2, 3}
Ceftriaxone	≤0.25	S ^{2, 3, 4}
Ciprofloxacin	≤0.25	S ^{2, 3}
Ertapenem	0.25	S ^{2, 3}
Gentamicin	≤1	S ^{2, 3}
Meropenem	≤0.25	S ^{2, 3}
Piperacillin-Tazobactam	≤4	S ^{2, 3}
Tobramycin	≤1	S ^{2, 3}

NOTES

AIM ² : Check culture for multiple morphologies due to possibility of a second species being present.
AIM ¹ : Identification results that are discordant with Gram stain should be confirmed with alternate method.
AIM ³ : Culture to agar media recommended to ensure the morphologic consistency of the ACCELERATE PHENO system identification
AIM ⁴ : For Enterobacterales, ceftriaxone SIR classification can be used to predict cefotaxime classification.



PennState Health



PennState
College of Medicine

How will the clinician view the result

Incorporate into the full susceptibility result

- Rapid AST results are interfaced directly to LIS
- Two AST reports
 - Rapid AST
 - Conventional AST
 - May not be needed if rapid AST has adequate antibiotics
- Hybrid AST report that contains both rapid and conventional AST results

Susceptibility

Hypothetical Hybrid Report

	Klebsiella pneumoniae +ANTIBIOTIC MIC		
Amoxicillin/CA	<=2 mcg/mL	Sensitive	
Ampicillin	16 mcg/mL	Resistant	
Aztreonam	<=1 mcg/mL	Sensitive	Rapid AST
Cefepime	<=1 mcg/mL	Sensitive	Rapid AST
Ceftriaxone	<=1 mcg/mL	Sensitive	Rapid AST
Ciprofloxacin	<=0.25 mcg/mL	Sensitive	
Ertapenem	<=0.5 mcg/mL	Sensitive	Rapid AST
Gentamicin	<=1 mcg/mL	Sensitive	Rapid AST
Levofloxacin	<=0.12 mcg/mL	Sensitive	
Meropenem	<=0.25 mcg/mL	Sensitive	Rapid AST
Piperacillin/tazobactam	<=4 mcg/mL	Sensitive	Rapid AST
Tetracycline	<=1 mcg/mL	Sensitive	
Trimethoprim/Sulfa	<=20 mcg/mL	Sensitive	



PennState Health



PennState
College of Medicine

What type of stewardship intervention is needed

- **Active**

- Laboratory phone call/ text message/ page/ email to a stewardship member on-call
 - i.e. Secure text to pharmacy resident between 7am-11pm 7 days a week
- EMR alert to stewardship or providers' inbox

- **Passive**

- Reporting comments

Klebsiella pneumoniae

This organism is carbapenem resistant. Infectious disease consult is recommended.

Pseudomonas aeruginosa

This organism is carbapenem resistant and cefepime susceptible. Infectious disease consult is recommended.



PennState Health



PennState
College of Medicine

What if the rapid and the conventional AST reports are discrepant

Rapid AST performance is good

- Categorical agreement
 - 85% - 95% compared to standard of care AST
- Essential agreement
 - 92% – 95% compared to standard of care AST
- Discordant analysis with reference BMD
 - ~50% of discordant results will agree with rapid AST

Comfort level of clinicians acting on results

- Performance of the rapid AST to your standard of care AST method
- Severity of illness of the patient
- Number of co-morbidities in a patient
- Additional pending laboratory tests and diagnostics

Tibetts et al, JCM 2022
Snyder et al, JCM 2024
Esse et al, JCM 2023
Unpublished, Penn State



PennState Health



PennState
College of Medicine

Additional Points to Consider

- Review your local antibiogram
 - Is there a high level of resistance in the population you serve
 - How often does the escalation of antibiotics occur
- Alignment with current institutional and professional treatment guidelines
 - 2024 IDSA guidance on the treatment of resistant gram-negative infections emphasizes the detection of specific resistance markers to guide treatment decisions
 - ESBL /CTX-M
 - KPC/ NDM/ OXA-48



PennState Health



PennState
College of Medicine

Thank You



PennState Health



PennState
College of Medicine