





**nmec**

RAMAN AND PL FOR NANOSCALE MATERIALS  
CHARACTERIZATION AND METROLOGY


# BASICS OF RAMAN SPECTROSCOPY

- Nonlinear optical technique  
→ measures light emitted with a (slightly) lower/higher energy than incident




- Si spectrum:
  - Rayleigh (linear scattering)
  - Stokes (emission of 1 phonon)
  - Anti-stokes (absorption of 1 phonon)

Raman spectrum of crystalline Si




# VALUES OF RAMAN SPECTROSCOPY

## Average composition




Rouchon J. Crys. Growth **392** 66 2014


## Average stress




## Geometry/CD



## Doping




## Different probing depths possible (profiling)




- Optical:
  - Fast
  - Non-contact
  - Good lateral resolution ( $\sim 1 \mu\text{m}$ )
  - Capable of probing in depth
- Sensitive to:
  - Material
  - Composition
  - Stress
  - Geometry

## Constant rise in measurement requests




# 30 YEARS OF RAMAN EXPERTISE AT IMEC

## Average stress measurements





De Wolf J, Appl. Phys. 71 898 (1992)  
De Wolf, Semicond. Sci. Technol. 11 139 (1996) – 720 citations

## Ultralow Si content





De Wolf ECS Trans. 86 397 (2018)

## Automated data treatment, fitting and plotting






## Stress/composition in 3D confined volumes




# STRESS AND CD ANALYSIS OF CONFINED VOLUMES

# NANOFOCUSING OF LIGHT

## ENABLING DEEP SUB-WAVELENGTH STRUCTURES TO STAND OUT




- Raman response dictated by geometric properties: critical dimension, pitch
- Allows for isolation of the signal of the ROI
- All strengths of Raman are re-enabled at the nanometer scale
- Opens path for Raman dimensional metrology







Kuhn FCMN 2013  
 Nuytten Appl. Phys. Lett. **106** 033107 (2015)  
 Bogdanowicz Appl. Phys. Lett. **108** 083106 (2016)  
 Nuytten APL Mat. **6** 058501 (2018)  
 Gawlik Appl. Phys. Lett. **113** 063103 (2018)

# CAVITY DEPTH




Etch time

## Nanosheet/Forksheet/CFET technology




CD~120 nm  
Pitch=520 nm




Cavity : 0 nm  
 $V_{SiGe} = V_0$

normalized integrated counts




--- V/V0  
—● Raman



Cavity: 28.7 nm  
 $V_{SiGe}/V_0 = 0.54$


## NANOSHEETS FORKSHEETS

CD~50 nm  
Pitch=90 nm




Cavity : 0 nm  
 $V_{SiGe} = V_0$

Normalized integrated counts




--- V/V0  
—● Raman

imec patent





Cavity ~ 20 nm  
 $V_{SiGe} \sim 0$

CD~120 nm  
Pitch=520 nm




Cavity : 0 nm  
 $V_{SiGe} = V_0$





Cavity: 28.7 nm  
 $V_{SiGe}/V_0 = 0.54$

## NANOSHEETS FORKSHEETS

CD~50 nm  
Pitch=90 nm




Cavity : 0 nm  
 $V_{SiGe} = V_0$




Cavity ~ 20 nm  
 $V_{SiGe} \sim 0$

# APPLICATION ON FULLY AUTOMATED TOOLS


## Cavity depth





## SiGe composition & strain








Strain in SiGe layer (x100%)



Strain in the Si sheet (x100%)



# UNIVERSAL APPLICATION SUITE




RECENT RAMAN DEVELOPMENTS:  
OFF-AXIS RAMAN  
AND  
USE OF LINEARIZED RADIAL POLARIZATION


# NANOSCALE STRESS MEASUREMENTS

## IMPROVE PRACTICALITY OR PRECISION OF ANISOTROPIC STRESS MEASUREMENTS


- Decoupling of  $\sigma_{11}$  and  $\sigma_{22}$  requires deconvolution of LO and TO
- Selection rules dictate TO only excited by  $E_z$
- High precision measurements = large out-of-plane polarization component  $E_z$
- Standard solution: oil immersion – no in-line potential




**Linear**



**Radial**






**Linearised Radial**



↕ Polarisation direction

# STRESS MEASUREMENT AND VALIDATION




Setup	TO/LO
Linear setup	$0.6 \pm 0.2$
Linearised radial setup	$1.2 \pm 0.3$



Linear	Linearised radial
$\sigma_x = -1.47 \pm 0.05 \text{ GPa}$	$\sigma_x = -1.45 \pm 0.03 \text{ GPa}$
$\sigma_y = -1.48 \pm 0.05 \text{ GPa}$	$\sigma_y = -1.45 \pm 0.03 \text{ GPa}$

**Biaxial in-plane stress**






Linear	Linearised radial
$\sigma_x = -0.45 \pm 0.05 \text{ GPa}$	$\sigma_x = -0.48 \pm 0.02 \text{ GPa}$
$\sigma_y = -2.29 \pm 0.05 \text{ GPa}$	$\sigma_y = -2.26 \pm 0.02 \text{ GPa}$




**Uniaxial stress**

**Improved precision in stress measurement**

# STRESS MEASUREMENTS WITH FCOR





- LO- and TO-active spectra are acquired independently
- In-plane stress extracted with similar precision
- Validated by OIR and NBD
- Truly non-invasive stress characterization



# TERS/TEPL FOR NANOSCALE RESOLUTION

# TIP ENHANCED RAMAN SPECTROSCOPY




## Concept


- Raman intensity  $\propto E^4$
- Incident laser excites a localized surface plasmon at the tip apex
- Plasmonically enhanced field increases the Raman intensity near the tip apex
- Intensity enhancement depends on tip properties, laser polarization and choice of substrate

## Spatial resolution


- Far field  $\sim \frac{0.61 \lambda}{N.A.}$
- Near field  $\sim$  Radius of tip apex




# TERS AFM topography




## Spectra




## Near field Intensity



—●— NF peak position    —■— fem



Equivalent stress (MPa)



# BACK-GATED CRYO-PHOTOLUMINESCENCE OF 2D MATERIALS

## NON-DESTRUCTIVE DEFECT CHARACTERIZATION?


Defects (e.g. S-vacancies) create shallow levels in the BG that can be studied through photoluminescence. For this we need low temperature (77K), and control over all other parameters affecting the PL:

Temperature


Excitation power

Back-gate bias / doping

Extrinsic & strain factors



Sample




~ 1cm

Laser spot location



1L MoS2/SiO2 at 77K



Built setup for combined back-gating and low-temperature PL, first electrostatic control and activation of defect signal on **imec synthetic material** has been demonstrated. Quantitative interpretation of results and link with electrical measurements needs development and data cycles.

### Two setups available at imec:

N2 sample holder  
(used on Raman tools)



Lifetime meter with cryo  
(less suited)



Steps  
to  
take

- Quantitative defect characterization from observed behavior
- Link with mobility (device processing)
- Establish connection with S-vacancies (HRTEM)

# ACKNOWLEDGEMENT



J. Bogdanowicz, Y. Oniki, K. Kenis, A. Minj, Z. Khan, S. Sergeant, V. Prabhakara, A. Franquet, V. Spampinato, T. Conard, I. Hofliijk, D. Van den Heuvel, H. Mertens, G. Eneman, L. Witters, S. Achra, S. Brems, C. Fleischmann, A.-L. Charley, P. van der Heide, P. Leray, I. De Wolf



J. Hung, R. Koret, D. Fishman, Y. Katz, N. Meir



umec

embracing a better life