




# Atom probe tomography using an extreme ultraviolet pulsed light source

Luis Miaja Avila, Benjamin Caplins, Jacob Garcia, Norman Sanford, Ann Chiamonti  
National Institute of Standards and Technology, Boulder, CO, USA

# NIST




Improve atom probe tomography instrumentation, experimental methods, and understanding of apex electrostatics.

**Thrust 1:** Build a new atom probe instrument optimized for the specific materials of interest to the semiconductor industry.


**Thrust 2:** Optimize atom probe measurements on commercial tools using specimens provided by the semiconductor industry.

**Thrust 3:** Develop modeling tools to simulate apex electrostatics on real-world specimens of interest to the semiconductor industry.




Atom probe microscopes operate in the ultraviolet region of the electromagnetic spectrum. The goal of Thrust 1 is to determine the optimal laser wavelength for the specific elements, materials, and structures of interest to the semiconductor industry.

## Atom Probe Tomography




Conventional APT measurements are performed with voltage pulses or laser pulses. Typically harmonics of a Nd-based laser are used having wavelengths around 532 nm, 355 nm, 257 nm (2.3 eV, 3.5 eV, 4.8 eV). For samples containing materials with large optical band gaps these photons are not efficiently absorbed and APT measurements can be challenging.

## Design of an EUV Atom Probe Microscope




## Characterization of a EUV Source for APT


### Characterizing the spectrum



EUV light requires a careful selection of optics to maximize beamline throughput. We use grazing incidence reflective optics to steer and focus the light which permit broadband reflection. We use ultrathin aluminum foils to filter the NIR light.




### Measuring the beam size and shape



To calculate the energy or photon fluence at the sample, the spot size at the sample plane must be known.

## EUV APT Data


### InGaN





QW (at. %)  
Ga = 40.5 ± 1.4  
N = 51.8 ± 1.4  
In = 7.7 ± 0.5

barrier (at. %)  
Ga = 48.1 ± 0.6  
N = 51.9 ± 0.6  
In = 0.0


buffer (at. %)  
Ga = 47.7 ± 0.4  
N = 50.1 ± 0.4  
In = 2.2 ± 0.1



InGaN quantum well sample, mass spectra, and EUV APT reconstruction. Using EUV light we obtained the correct indium apparent composition values.



### GaN



GaN EUV APT mass spectra. Apparent compositional measurements show clear biases in gallium and nitrogen compositions, similar to NUV APT. These measurements suggest the electric field at the tip is the main cause of compositional bias in GaN.

### References:

Caplins et al., Atom Probe Tomography Using an Extreme Ultraviolet Trigger Pulse. Rev. Sci. Instrum. 94, 093704 (2023) DOI: 10.1063/5.0160797  
Garcia et al., A Comprehensive Examination of Al<sub>2</sub>O<sub>3</sub> Using EUV and NUV LAPT. M&M 29, 83 (2023) DOI: 10.1093/micmic/ozad067.033  
Miaja et al., EUV Radiation Pulsed APT of III-Nitride Semiconductor Materials. JPCC 125, 2626 (2021) DOI: 10.1021/acs.jpcc.0c08753  
Chiamonti et al. Field Ion Emission in an APT Triggered by Coherent EUV Light. M&M 26, 258 (2020) DOI: 10.1017/S1431927620000203

We thank Daniel Hickstein, Seth Cousin, and Clayton Bargsten for their help in the design of the EUV beamline. We gratefully acknowledge financial support from the NIST CHIPS Program and NIST Innovations in Measurement Science Program and technical support from CAMECA Instruments through a cooperative research and development agreement.